

Datasheet

GPR115 polyclonal antibody

Catalog Number: PAB26314

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GPR115.

Immunogen: A synthetic peptide corresponding to 16 amino acids at C-terminus of human GPR115.

Host: Rabbit

Reactivity: Bovine, Gorilla, Horse, Human, Monkey, Mouse

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: BLAST analysis of the peptide immunogen showed no homology with other human proteins, except GPR113 (69%).

Form: Liquid

Purification: Immunoaffinity chromatography

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (4-10 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 221393

Gene Symbol: GPR115

Gene Alias: FLJ38076, PGR18