



Certificate of Conformity

VWR CELL CULTURE FLASKS, STANDARD LINE

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Cap	Capacity(ml)	Culture area(cm ²)	Packed(/bag)	Pk
734-2309	Vented	25	12.5	10	200
734-2310	Plug seal	25	12.5	10	200
734-2311	Vented	50	25	10	200
734-2312	Plug seal	50	25	10	200
734-2313	Vented	250	75	5	100
734-2314	Plug seal	250	75	5	100
734-2315	Vented	600	182.5	5	40
734-2316	Plug seal	600	182.5	5	40
734-2600	Vented	850	300	3	18
734-2601	Plug seal	850	300	3	18

Material:

Flask: Polystyrene PS1050 Color: Clear

Cap(734-2309 to 734-2316): High Density Polyethylene 2911

Color: Purple 259

Cap(734-2600 & 734-2601): High Density Polyethylene HMA025

Color: Purple 259

Membrane: 0.22µm Hydrophobic PP

Features:

- Temperature range: -20°C to +50°C
- Autoclavability: No
- Non-pyrogenic
- Products does not contain material or component containing Phtalane, Biocide, Oleamide, Erucamide, Stearamide, Silicone, Melamine, Latex and TSE/BSE
- Shelf life: 3 years after month of production


Sterility:

- Sterilized by gamma irradiation SAL10⁻⁶




Certificate of Conformity

734-2309 & 734-2310




734-2311 & 734-2312








Certificate of Conformity


734-2313 & 734-2314





734-2315 & 734-2316




vented cap




plug seal cap



vented cap




plug seal cap





Certificate of Conformity

734-2600 & 734-2601



Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2016**