



Certificate of Conformity

VWR CELL TISSUE CULTURE FLASK, TREATED

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Description	Working Volume	Growth Area	Qty /Pk	Qty/ Case
734-2810	600 ml, Treated for Increased Cell Attachment, w/plug seal cap, Sterilized	400 ml	182 cm ²	20	40
734-2811	600 ml, Treated for Increased Cell Attachment, w/vent cap, Sterilized	400 ml	182 cm ²	20	40
734-2812	850 ml, Treated for Increased Cell Attachment, w/plug seal cap, Sterilized	800 ml	300 cm ²	9	18
734-2813	850 ml, Treated for Increased Cell Attachment, w/vent cap, Sterilized	800 ml	300 cm ²	9	18

Materials:

Flask: GPPS (General polystyrene) Colour: Clear

Cap: HDPE(High Density Polyethylene) Colour: Purple 259

Membrane: 0.22µm Hydrophobic PP (Only use for Vented Cap)




Certificate of Conformity

734-2810 and 734-2811



3: Dimensions (Unit: mm)






Certificate of Conformity

734-2812 and 734-2813



Dimensions (Unit: mm)



Certificate of Conformity

Other Characteristics:

- Temperature: -20°C~+50°C
- Sterilized by gamma irradiation SAL10⁻⁶(ISO11137)
- Non-pyrogenic
- Bulk packed in zip-Lock bag
- Autoclavability: No
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Product is latex free.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2017**